



MICROMATTER Technologies Inc.
8333-130th Street, Unit#1, Surrey BC
V3W 7X4 Canada
www.micromatter.com

Product Information and Frequently Asked Questions

Special XRF Calibration Materials

Dual Element and Variable Stoichiometry Standards / Stoichiometric factors for dual element standards

April 11, 2018

Micromatter has manufactured calibration and reference materials for several decades. Our products contain highest quality materials, such as ultra-pure metals or stable inorganic compounds.

What are dual element standards?

Micromatter dual element standards contain compounds of well-defined composition. The area weight noted on the standard and the certification sheet refers to the entire compound. The cation and anion contents can be calculated from the molecular mass of the compound using the stoichiometric ratios of the two components.

Dual element standards are an economic alternative to single element standards as they can be used for calibrating an XRF spectrometer in different energy regions. For example, zinc telluride (ZnTe) can be used to calibrate both the zinc and the tellurium peak. It is not necessary to purchase a separate pure tellurium standard.

Why are some Micromatter standards marked as having 'variable stoichiometry'?

Micromatter produces a number of dual element standards, such as gallium phosphide or gallium arsenide, by evaporation of solid binary compounds onto the desired backing material (Nuclepore membrane or Mylar). Our manufacturing technique follows an established process commonly used in the semiconductor industry.

It has been evaluated that the resulting deposits have a different composition than the starting material. This composition is variable and depends on various factors during the material deposition. Micromatter has recently verified the accuracy of these particular standards by independent chemical analysis. Therefore, any standards that are certified to have a variable stoichiometry are accompanied with individual concentrations/loadings specified for both of the elements in the compound. It is crucial that for compounds with variable stoichiometry, the loadings specified by Micromatter be used, rather than those calculated using stoichiometric factors.

Stoichiometric Factors for dual-element standards with non-variable stoichiometry

The following table lists convenient conversion factors for the calculation of elemental weights from given area weights.

Note: The stoichiometric factor for all pure element standards is 1.00.

Example:

What is the area weight of scandium in my ScF₃ standard, if its certified weight is 45.8 μg/cm²?

$$45.8 \mu\text{g}/\text{cm}^2 \times 0.441 \text{ (tabulated factor for Sc)} = 20.16 \mu\text{g}/\text{cm}^2$$

Element	Compound	Factor		Element	Compound	Factor
Li	LiF	0.250		La	LaF ₃	0.709
Na	NaCl	0.393		Ce	CeF ₃	0.711
Mg	MgF ₂	0.390		Pr	PrF ₃	0.712
Si	SiO	0.637		Nd	NdF ₃	0.717
P	GaP	0.308		Sm	SmF ₃	0.725
S	CuS _x	*		Eu	EuF ₃	0.727
Cl	KCl	0.476		Gd	GdF ₃	0.734
	NaCl	0.607		Tb	TbF ₃	0.736
K	KCl	0.524		Dy	DyF ₃	0.740
	KI	0.236		Ho	HoF ₃	0.743
Ca	CaF ₂	0.513		Er	ErF ₃	0.746
Sc	ScF ₃	0.441		Tm	TmF ₃	0.748
Zn	ZnTe	0.339		Yb	YbF ₃	0.752
Ga	GaP	0.692		Lu	LuF ₃	0.754
	GaAs	0.482				
As	GaAs	0.518		W	WO ₃	0.793
Br	CsBr	0.375		Hg	AgHg	*
Rb	RbI	0.402		Tl	TlCl	0.852
Sr	SrF ₂	0.698				
Y	YF ₃	0.609				
Zr	ZrF ₄	0.546				
Nb	Nb ₂ O ₃	0.795				
Mo	MoO ₃	0.667				
Cd	CdSe	0.587				
I	RbI	0.598				
	KI	0.764				
Cs	CsBr	0.625				
Ba	BaF ₂	0.783				

*Composition varies. See Certification Sheet for details.